

Gastrin 17 (G-17) Mouse mAb

Cat. No: V4001
Name: Gastrin 17 (G-17) Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No:
Host/Source: Mouse
Subtype: IgG
Reactivity: (predicted: Human)
Format: Liquid
Concentration: $\geq 1\text{mg/ml}$
Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications: CLIA
Antigen: Recombinant human Gastrin 17
Pairing Information:

CL		Detection	
		V4001	V4002
Capture	V4001	-	+++
	V4002	-	-