

## Gastrin 17 (G-17) Mouse mAb

Cat. No: V4002  
Name: Gastrin 17 (G-17) Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No:  
Host/Source: Mouse  
Subtype: IgG  
Reactivity: (predicted:Human)  
Format: Liquid  
Concentration:  $\geq 1$ mg/ml  
Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.  
Applications: CLIA  
Antigen: Recombinant human Gastrin 17  
Pairing Information:

CL		Detection	
		V4001	V4002
Capture	V4001	-	+++
	V4002	-	-