

Human IgE Mouse mAb

Cat. No: V5402
Name: Human IgE Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5E9
Host/Source: Mouse
Subtype: IgG2b
Reactivity: Human
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Purified native protein ([bs-2226P](#))
Pairing Information:

ELISA		Detection		
		V5402	V5404	V5406
Capture	V5402	-	-	-
	V5404	+++	-	-
	V5406	+++	-	-
ICA		Detection		
		V5402	V5404	V5406
Capture	V5402	-	-	-
	V5404	+++	-	-
	V5406	+++	-	-