
Creatine Kinase MB (CK-MB) Mouse mAb

Cat. No:	V1002
Name:	Creatine Kinase MB (CK-MB) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7C9
Host/Source:	Mouse
Subtype:	IgG1
Format:	Liquid
Concentration:	≥0.6 mg/ml
Buffer:	potassium phosphate buffer and sodium chloride. No preservatives are added.
Applications:	ELISA, ICA
Antigen:	Recombinant human Creatine Kinase MB (CK-MB) (bs-41114P)