

bs-6081R

[Primary Antibody]

PDGFRL Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 5157 Target: PDGFRL Immunogen: KLH conjugated synthetic peptide derived from human PDGFRL: 61-150/375. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: This protein is expressed in colon, lung and liver. Defects in PDGFRL are associated with colorectal cancer.	Isotype: IgG SWISS: Q15198	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Dog, Horse) Predicted MW.: 40 kDa Subcellular Location: Secreted
--	--	---