bs-19139R

[Primary Antibody]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

MYBPHL Rabbit pAb

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 343263 SWISS: A2RUH7

Target: MYBPHL

Immunogen: KLH conjugated synthetic peptide derived from human MYBPHL:

231-330/354.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit)

Predicted MW.: 39 kDa

Subcellular Location: Cytoplasm