

bs-14124R

[Primary Antibody]

EZH1P Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Human)
GeneID: 340602	SWISS: Q86X51	
Target: EZHIP		
Immunogen: KLH conjugated synthetic peptide derived from human CXorf67: 301-400/503.		Predicted MW.: 52 kDa
Purification: affinity purified by Protein A		Subcellular Location: Nucleus
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The function of this protein is not yet known.		