

bs-15298R

[Primary Antibody]

C8ORF59 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 401466	SWISS: Q8N0T1	Reactivity: (predicted: Human, Mouse, Rat, Rabbit)
Target: C8ORF59		
Immunogen: KLH conjugated synthetic peptide derived from human C8ORF59: 1-100/100.		
Purification: affinity purified by Protein A		Predicted MW.: 11 kDa
Concentration: 1mg/ml		Subcellular Location: Cytoplasm ,Nucleus
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=6.1]** Pan Hongyu. et al. RBIS regulates ribosome biogenesis to affect progression in lung adenocarcinoma. J TRANSL MED. 2024 Dec;22(1):1-17 IF ;Human. 39722009