

**bs-2032R****[ Primary Antibody ]****PMP2 Rabbit pAb****BioSS**  
**ANTIBODIES**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**DATASHEET****Host:** Rabbit**Isotype:** IgG**Clonality:** Polyclonal**GeneID:** 5375**Target:** PMP2**Immunogen:** KLH conjugated synthetic peptide derived from human PMP2: 31-132/132.**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


**Background:** PMP2 is a small basic protein, belonging to the fatty-acid binding protein (FABP) family and is one of the major proteins of peripheral myelin. It appears to be involved in the transport of fatty acids or the metabolism of myelin lipids, and has been shown to have lipid-binding activity. PMP2 protein may have an important role in the organization of compact myelin.**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)


IF (1:100-500)

Flow-Cyt (1ug/Test)


ICC/IF (1:100)

**Reactivity:** Human, Mouse, Rat  
(predicted: Rabbit, Pig, Cow)**Predicted MW.:** 15 kDa**Subcellular Location:** Cytoplasm**VALIDATION IMAGES**


Tissue/cell: Mouse embryos; 4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01M, pH 6.0), Boiling bathing for 15min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30min; Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min; Incubation: Anti-phospho-KLF5(Ser275) Polyclonal Antibody, Unconjugated(bs-2020R) 1:200, overnight at 4°C, followed by conjugation to the secondary antibody(SP-0023) and DAB(C-0010) staining



SH-SY5Y cell; 4% Paraformaldehyde-fixed; Triton X-100 at room temperature for 20 min; Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min; Antibody incubation with (PMP2) polyclonal Antibody, Unconjugated (bs-2032R) 1:100, 90 minutes at 37°C; followed by a conjugated Goat Anti-Rabbit IgG antibody at 37°C for 90 minutes, DAPI (blue, C02-04002) was used to stain the cell nuclei.



Tissue/cell: Rat brain; 4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01M, pH 6.0), Boiling bathing for 15min; Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min; Incubation: Anti-PMP2 Polyclonal Antibody, Unconjugated(bs-2032R) 1:200, overnight at 4°C; The secondary antibody was Goat Anti-Rabbit IgG, Cy3 conjugated(bs-0295G-Cy3) used at 1:200 dilution for 40 minutes at 37°C. DAPI(5ug/ml, blue, C-0033) was used to stain the cell nuclei



Blank control (black line) :SH-SY5Y. Primary Antibody (green line): Rabbit Anti-PMP2 antibody (bs-2032R) Dilution:1ug/Test; Secondary Antibody (white blue line) : Goat anti-rabbit IgG-AF488 Dilution: 0.5ug/Test. Isotype control (orange line) : Normal Rabbit IgG Protocol The cells were fixed with 4% PFA

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

(10min at room temperature)and then  
permeabilized with 90% ice-cold methanol for  
20 min at -20°C, The cells were then incubated in  
5%BSA to block non-specific protein-protein  
interactions for 30 min at room temperature  
.Cells stained with Primary Antibody for 30 min  
at room temperature. The secondary antibody  
used for 40 min at room temperature.  
Acquisition of 20,000 events was performed.