

bs-13552R

[**Primary Antibody**]

Bioss[™]
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

ZBED2 Rabbit pAb

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 79413</p> <p>Target: ZBED2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZBED2: 1-100/218.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ZBED3 is an Axin-binding protein involved in Wnt/beta-catenin signaling modulation. ZBED4 expression has been shown in human and mouse retinas where it is thought to act as a regulatory protein in cone photoreceptors and Müller cells. The functions of ZBED2 and ZBED5 have yet to be elucidated.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9BTP6</p>	<p>Applications: ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Sheep, Cow)</p> <p>Predicted MW.: 25 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--	--