

bsm-2079M

[Primary Antibody]

Cecropin-A Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG

Clonality: Monoclonal

CloneNo.: 1D8

Target: Cecropin-A

Immunogen: KLH conjugated synthetic peptide derived from Hyalophora cecropia: 1-37/41.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: Size : 50ul/100ul/200ul
0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Size : 200ug (PBS only)
0.01M PBS
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Hyalophora cecropia)

**Predicted
MW.:** 4 kDa