
AHD Mouse mAb

Catalog Number: bsm-20932M

Target Protein: AHD

Concentration: 1mg/ml

Form: Liquid

Host: Mouse

Clonality: Monoclonal

Clone No.: 5F8

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:AHD)

Purification: affinity purified by Protein G

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.