

bsm-52960R

[Primary Antibody]

CRISPR-Cas9 Recombinant Rabbit mAb

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:50-200) IHC-F (1:50-200) IF (1:50-200) Reactivity: Streptococcus pyogenes Predicted MW.: 160 kDa
Clonality: Recombinant	CloneNo.: 6C1	
Target: CRISPR-Cas9		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.02% Proclin300. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— VALIDATION IMAGES —



Sample: Lane 1: Recombinant CRISPR-Cas9 protein, C-His(bs-41498P) Primary: Anti-CRISPR-Cas9 (bsm-52960R) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 160 kDa
Observed band size: 182 kDa