
SEM Mouse mAb

Catalog Number: bsm-20934M

Target Protein: SEM

Concentration: 1mg/ml

Form: Liquid

Host: Mouse

Clonality: Monoclonal

Clone No.: 3H11

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:SEM)

Purification: affinity purified by Protein G

Storage: 0.01M TBS (pH7.4).

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.