

**bsm-20934M**

**[ Primary Antibody ]**

## SEM Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG

**Clonality:** Monoclonal

**CloneNo.:** 3H11

**Target:** SEM

**Purification:** affinity purified by Protein G

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4).  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated  
freeze/thaw cycles.

**Applications:** ELISA (1:5000-10000)

**Reactivity:** (predicted: SEM)