bsm-20934M

## [ Primary Antibody ]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## Nitrofurazone Mouse mAb

- DATASHEET -

Host: Mouse Isotype: IgG Clonality: Monoclonal CloneNo.: 3H11

Target: Nitrofurazone

**Purification:** affinity purified by Protein G

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4). Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: ELISA (1:5000-10000)

Reactivity: Species independent