

bsm-20934M

[Primary Antibody]

SEM Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Clonality: Monoclonal

Target: SEM

Purification: affinity purified by Protein G

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4).
Shipped at 4°C. Store at -20°C for one year. Avoid repeated
freeze/thaw cycles.

Isotype: IgG

CloneNo.: 3H11

Applications: ELISA (1:5000-10000)

Reactivity: Species independent