

RBP4 Mouse mAb

Cat. No: V6202
Name: RBP4 Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 7D3
Host/Source: Mouse
Subtype: IgG2b
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human RBP4 Protein ([bs-43011P](#))
Pairing Information:

ELISA		Detection				
		V6201	V6202	V6203	V6204	V6205
Capture	V6201	-	-	+	-	-
	V6202	-	-	-	-	-
	V6203	-	+++	-	-	-
	V6204	-	-	+++	-	-
	V6205	-	+	-	-	-
ICA		Detection				
		V6201	V6202	V6203	V6204	V6205
Capture	V6201	-	+++	+	-	-
	V6202	-	-	-	-	-
	V6203	+++	-	-	-	-
	V6204	-	-	-	-	-
	V6205	-	-	-	-	-