

Cystatin-C Mouse mAb

Cat. No:	V0506
Name:	Cystatin-C Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C1
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human Cys-C Protein (bs-41175P)
Pairing Information:	

ELISA		Detection			
		V0503	V0505	V0506	V0507
Capture	V0503	-	-	-	-
	V0505	-	-	-	-
	V0506	+++	-	-	-
	V0507	+++	-	-	-
ICA		Detection			
		V0503	V0505	V0506	V0507
Capture	V0503	-	-	+++	-
	V0505	-	-	+++	-
	V0506	-	-	-	-
	V0507	-	-	-	-