

---

## Human S100B Mouse mAb

|               |   |
|---------------|---|
| Cat. No:      | V2506   |
| Name:         | Human S100B Mouse mAb   |
| Storage:      | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:   | Monoclonal  |
| Clone No:     | 3B8   |
| Host/Source:  | Mouse   |
| Subtype:      | IgG1  |
| Reactivity:   | (predicted:&nbsp;Human)   |
| Format:       | Liquid  |
| Buffer:       | 0.01M PBS (pH7.4).  |
| Applications: | ELISA   |
| Antigen:      | Recombinant Human S100B Protein ( <a href="#">bs-41191P</a> )                   |