

Human S100B Mouse mAb

Cat. No:	V2504
Name:	Human S100B Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5A9
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant human S100B protein (bs-41191P)
Pairing Information:	

ELISA		Detection			
		V2501	V2503	V2504	V2505
Capture	V2501	-	-	-	-
	V2503	+++	-	-	+
	V2504	-	+	-	-
	V2505	+	+	-	-
ICA		Detection			
		V2501	V2503	V2504	V2505
Capture	V2501	-	-	-	-
	V2503	+	-	+++	-
	V2504	-	-	-	-
	V2505	-	-	-	-