
COT Mouse mAb

Catalog Number: bsm-3632M

Target Protein: COT

Concentration: 1mg/ml

Form: Liquid

Host: Mouse

Clonality: Monoclonal

Clone No.: 7G1

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:COT)

Purification: affinity purified by Protein G

Storage: 0.01M TBS (pH7.4).

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.