

bsm-4474M

[Primary Antibody]

BZO Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG

Clonality: Monoclonal

CloneNo.: 12D6

Target: BZO

Purification: affinity purified by Protein G

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4).
Shipped at 4°C. Store at -20°C for one year. Avoid repeated
freeze/thaw cycles.

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: BZO)