

**bsm-4476M**

**[ Primary Antibody ]**

## MDMA Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Mouse	<b>Isotype:</b> IgG	<b>Applications:</b> ELISA (1:5000-10000) <b>Reactivity:</b> Species independent
<b>Clonality:</b> Monoclonal	<b>CloneNo.:</b> 5F2	
<b>Target:</b> MDMA		
<b>Purification:</b> affinity purified by Protein G		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4). Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		