

Prostate Specific Antigen Mouse mAb

Cat. No: V1704
Name: Prostate Specific Antigen Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5F2
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Prostate Specific Antigen (PSA) Antigen, Native ([bs-41609P](#))
Pairing Information:

ELISA		Detection	
		V1703	V1704
Capture	V1703	-	+++
	V1704	-	-
ICA		Detection	
		V1703	V1704
Capture	V1703	-	-
	V1704	+++	-