

---

## Interleukin-1 receptor-like 1/ST2 Mouse mAb

Cat. No:	V3701
Name:	Interleukin-1 receptor-like 1/ST2 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1E3
Host/Source:	Mouse
Subtype:	IgG1
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human ST2 Protein ( <a href="#">bs-43005P</a> )