
Interleukin-1 receptor-like 1/ST2 Mouse mAb

| | |
|---------------|---|
| Cat. No: | V3703 |
| Name: | Interleukin-1 receptor-like 1/ST2 Mouse mAb |
| Storage: | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type: | Monoclonal |
| Clone No: | 2B2 |
| Host/Source: | Mouse |
| Subtype: | IgG1 |
| Format: | Liquid |
| Buffer: | 0.01M PBS (pH7.4). |
| Applications: | ELISA |
| Antigen: | Recombinant Human ST2 Protein (bs-43005P) |