
NRH2 Rabbit pAb

Catalog Number: bs-7672R

Target Protein: NRH2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Mouse, Rat, Pig, Sheep, Dog, Horse)

Predicted MW: 24 kDa

Source: KLH conjugated synthetic peptide derived from rat NRH2: 51-150/228.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: A single transmembrane protein highly expressed in embryonic tissues and involved in inducing apoptosis.