

**bs-7672R**

**[ Primary Antibody ]**

## NRH2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** NRH2

**Immunogen:** KLH conjugated synthetic peptide derived from rat NRH2: 51-150/228.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** A single transmembrane protein highly expressed in embryonic tissues and involved in inducing apoptosis.

**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

**Reactivity:** (predicted: Mouse, Rat, Pig, Sheep, Dog, Horse)

**Predicted MW.:** 24 kDa

**Subcellular Location:** Cell membrane ,Nucleus