

Azurocidin/HBP Mouse mAb

Cat. No: V4604
Name: Azurocidin/HBP Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 3F3
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human Azurocidin Protein ([bs-43022P](#))
Pairing Information:

ELISA		Detection				
		V4601	V4602	V4603	V4604	V4605
Capture	V4601	-	+	-	-	-
	V4602	-	-	-	-	-
	V4603	-	-	-	+++	-
	V4604	+++	-	-	-	-
	V4605	-	-	-	-	-
ICA		Detection				
		V4601	V4602	V4603	V4604	V4605
Capture	V4601	-	-	-	-	-
	V4602	-	-	-	-	-
	V4603	-	-	-	-	-
	V4604	-	-	+++	-	-
	V4605	-	-	-	+++	-