

## Azurocidin/HBP Mouse mAb

Cat. No: V4605  
Name: Azurocidin/HBP Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 6C7  
Host/Source: Mouse  
Subtype: IgG1  
Reactivity: (predicted:&nbsp;Human)  
Format: Liquid  
Buffer: 0.01M PBS (pH7.4).  
Applications: ICA  
Antigen: Recombinant Human Azurocidin Protein ([bs-43022P](#))  
Pairing Information:

| ELISA   |       | Detection |       |       |       |       |
|---------|-------|-----------|-------|-------|-------|-------|
|         |       | V4601     | V4602 | V4603 | V4604 | V4605 |
| Capture | V4601 | -         | +     | -     | -     | -     |
|         | V4602 | -         | -     | -     | -     | -     |
|         | V4603 | -         | -     | -     | +++   | -     |
|         | V4604 | +++       | -     | -     | -     | -     |
|         | V4605 | -         | -     | -     | -     | -     |
| ICA     |       | Detection |       |       |       |       |
|         |       | V4601     | V4602 | V4603 | V4604 | V4605 |
| Capture | V4601 | -         | -     | -     | -     | -     |
|         | V4602 | -         | -     | -     | -     | -     |
|         | V4603 | -         | -     | -     | -     | -     |
|         | V4604 | -         | -     | +++   | -     | -     |
|         | V4605 | -         | -     | -     | +++   | -     |