

---

## **beta-Amyloid(1-42) Mouse mAb**

Cat. No:	V5302
Name:	beta-Amyloid(1-42) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C8
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated synthetic peptide derived from human beta-Amyloid(1-42) ( <a href="#">bs-0107P</a> )
Pairing Information:	V5301&V5303-V5202