

## beta-Amyloid (1-40) Mouse mAb

Cat. No:	V5310
Name:	beta-Amyloid (1-40) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3D12
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted:&nbsp;Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	KLH conjugated synthetic peptide derived from human beta-Amyloid(1-40) ( <a href="#">bs-0106P</a> )
Pairing Information:	

ELISA		Detection								
		V5304	V5305	V5306	V5307	V5308	V5309	V5310	V5311	V5312
Capture	V5304	-	-	-	-	-	-	-	-	-
	V5305	+	-	-	-	-	-	-	-	-
	V5306	+	-	-	-	-	-	-	-	-
	V5307	+	-	-	-	-	-	-	-	-
	V5308	+	-	-	-	-	-	-	-	-
	V5309	-	+	-	-	-	-	-	+	-
	V5310	-	-	-	-	-	-	-	-	-
	V5311	-	-	-	-	-	-	-	-	-
	V5312	-	-	-	-	-	-	+	+	-