

Human Tau-4 Mouse mAb

Cat. No:	V5203
Name:	Human Tau-4 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2A6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human Tau-4 Protein (bs-41332P)
Pairing Information:	

ELISA		Detection																
		V5203	V5206	V5214	V5215	V5216	V5218	V5219	V5232	V5233	V5234	V5236	V5237	W5201	V5209	V5210	V5211	V5226
Capture	V5203	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	V5204	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-	-	-
	V5205	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	V5206	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	V5207	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	V5212	+	-	-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
	V5213	+	+	+	+	+	-	-	-	-	-	-	-	-	+	+	+	-
	V5215	+	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
	V5216	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	V5217	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	V5218	-	-	-	-	-	-	+	-	-	-	-	-	+	-	-	-	-
	V5220	-	-	-	-	-	-	-	-	-	+	-	-	+	-	-	-	-
	V5231	-	-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-
	V5233	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+++
	V5234	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
	V5235	-	-	-	-	-	-	+	+	-	+	-	-	-	-	-	-	-
	V5236	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	+++
	W5201	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
	V5225	-	+	+	-	-	+	-	-	-	-	+	-	-	-	-	-	-
	V5227	+	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
	V5222	-	-	+	-	-	+	-	-	-	-	-	+	-	-	-	-	-
	V5223	-	-	+	-	-	+	-	-	-	-	-	+	-	-	-	-	-

PRODUCT SPECIFIC PUBLICATIONS

[IF=5.7] Jihong Sun. et al. A Fully Automated Paper-Based Smartphone-assisted Microfluidic Chemiluminescence Sample-to-Result Immunoassay Platform. ANAL CHIM ACTA. 2025 Apr;;344013 ; . 40374237