

Human Tau Rabbit pAb

Cat. No: W5201
 Name: Human Tau Rabbit pAb
 Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
 Clone type: Polyclonal
 Clone No:
 Host/Source: Rabbit
 Subtype: IgG
 Reactivity: (predicted: Human)
 Purity: >90% as determined by SDS-PAGE
 Format: Liquid
 Buffer: 0.01M PBS (pH7.4).
 Applications: ELISA
 Antigen: Recombinant Human Tau-4 Protein ([bs-41332P](#))

Pairing Information:

ELISA	Detection																
	V5203	V5206	V5214	V5215	V5216	V5218	V5219	V5232	V5233	V5234	V5236	V5237	W5201	V5209	V5210	V5211	V5226
Capture	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-