
phospho-Tau (Thr217) Mouse mAb

Cat. No:	V5208
Name:	phospho-Tau (Thr217) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3C4
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr217