

NF-L Mouse mAb

Cat. No:	V5501
Name:	NF-L Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3H9
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Mouse, Rat (predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human NF-L Protein (bs-41217P)
Pairing Information:	

ELISA		Detection									
		V5501	V5502	V5503	V5504	V5505	V5506	V5507	V5508	V5509	V5510
Capture	V5501	-	-	-	-	+	-	-	-	-	-
	V5502	-	-	-	-	-	-	-	-	-	-
	V5503	-	-	-	-	+	-	-	-	-	-
	V5504	-	+	-	-	-	-	-	-	-	-
	V5505	-	-	-	+	-	+	-	-	-	-
	V5506	-	-	-	-	-	-	-	-	-	-
	V5507	-	+	-	-	-	-	-	-	-	-
	V5508	-	-	-	-	-	-	-	-	+	-
	V5509	-	+++	-	-	-	-	-	-	-	-
	V5510	-	+++	-	-	-	-	-	-	-	-