

bs-9293R

[Primary Antibody]

BNC2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:50-200) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Chicken, Dog, Horse) Predicted MW.: 122 kDa Subcellular Location: Cell membrane ,Cytoplasm ,Nucleus
Clonality: Polyclonal		
GeneID: 54796	SWISS: Q6ZN30	
Target: BNC2		
Immunogen: KLH conjugated synthetic peptide derived from human BNC2/Basonuclin 2: 601-700/1099.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Probable transcription factor specific for skin keratinocytes. May play a role in the differentiation of spermatozoa and oocytes. There are two isoforms.		

— SELECTED CITATIONS —

- **[IF=4.174]** Zhang, Zhao. et al. Development of a CAFs-related gene signature to predict survival and drug response in bladder cancer. Hum Cell. 2022 Jan;;1-16 IHC ;Human. 35044630