
NKX1-2 Rabbit pAb

Catalog Number: bs-24100R

Target Protein: NKX1-2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:5000-10000)

Reactivity: (predicted:Mouse)

Predicted MW: 32 kDa

Entrez Gene: 20231

Swiss Prot: P42580

Source: KLH conjugated synthetic peptide derived from mouse NKX1-2: 1-100/310.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: NKX1-2 (NK1 Homeobox 2) is a Protein Coding gene. GO annotations related to this gene include sequence-specific DNA binding. An important paralog of this gene is NKX1-1.