

**bs-24110R**

**[ Primary Antibody ]**

## PRR19 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                                                        |                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 623131</p> <p><b>Target:</b> PRR19</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from mouse PRR19: 81-180/366.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br/>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> The specific function of PRR19 remains unknown. There are 2 isoforms produced by alternative splicing.</p> | <p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> B2RW88</p> | <p><b>Applications:</b> <b>WB</b> (1:500-2000)<br/><b>IHC-P</b> (1:100-500)<br/><b>IHC-F</b> (1:100-500)<br/><b>IF</b> (1:100-500)<br/><b>ICC/IF</b> (1:100-500)</p> <p><b>Reactivity:</b> (predicted: Mouse, Rat)</p> <p><b>Predicted MW.:</b> 39 kDa</p> <p><b>Subcellular Location:</b> Nucleus</p> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|