

## AFP Mouse mAb

Cat. No: V2003  
 Name: AFP Mouse mAb  
 Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
 Clone type: Monoclonal  
 Clone No: 2D7  
 Host/Source: Mouse  
 Subtype: IgG1  
 Reactivity: (predicted:&nbsp;Human)  
 Format: Liquid  
 Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.  
 Applications: ELISA, ICA  
 Antigen: Recombinant Human AFP Protein ([bs-43140P](#))  
 Pairing Information:

ELISA		Detection				
		V2003	V2004	V2005	V2006	V2007
Capture	V2003	-	+++	-	-	-
	V2004	-	-	-	-	-
	V2005	-	+++	-	-	-
	V2006	-	-	-	-	-
	V2007	-	-	-	-	-
ICA		Detection				
		V2003	V2004	V2005	V2006	V2007
Capture	V2003	-	+++	-	-	-
	V2004	-	-	-	-	-
	V2005	-	-	-	-	-
	V2006	-	-	-	-	-
	V2007	-	-	-	+++	-