

AFP Mouse mAb

Cat. No: V2004
Name: AFP Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 1B1
Host/Source: Mouse
Subtype: IgG2a
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications: ELISA, ICA
Antigen: Recombinant Human AFP Protein ([bs-43140P](#))
Pairing Information:

ELISA		Detection				
		V2003	V2004	V2005	V2006	V2007
Capture	V2003	-	+++	-	-	-
	V2004	-	-	-	-	-
	V2005	-	+++	-	-	-
	V2006	-	-	-	-	-
	V2007	-	-	-	-	-
ICA		Detection				
		V2003	V2004	V2005	V2006	V2007
Capture	V2003	-	+++	-	-	-
	V2004	-	-	-	-	-
	V2005	-	-	-	-	-
	V2006	-	-	-	-	-
	V2007	-	-	-	+++	-

PRODUCT SPECIFIC PUBLICATIONS

[IF=34.9] Yi Yakun. et al. Nanopore-based enzyme-linked immunosorbent assay for cancer biomarker detection. NAT NANOTECHNOL. 2025 May;;1-8 ELISA ; Human . 40369343