

## AFP Mouse mAb

**Cat. No:** V2006  
**Name:** AFP Mouse mAb  
**Storage:** Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
**Clone type:** Monoclonal  
**Clone No:** 1B5  
**Host/Source:** Mouse  
**Subtype:** IgG1  
**Reactivity:** (predicted:&nbsp;Human)  
**Format:** Liquid  
**Buffer:** 0.01M PBS (pH7.4) with 0.02% Proclin300.  
**Applications:** ICA  
**Antigen:** Recombinant Human AFP Protein ([bs-43140P](#))  
**Pairing Information:**

ELISA		Detection				
		V2003	V2004	V2005	V2006	V2007
Capture	V2003	-	+++	-	-	-
	V2004	-	-	-	-	-
	V2005	-	+++	-	-	-
	V2006	-	-	-	-	-
	V2007	-	-	-	-	-
ICA		Detection				
		V2003	V2004	V2005	V2006	V2007
Capture	V2003	-	+++	-	-	-
	V2004	-	-	-	-	-
	V2005	-	-	-	-	-
	V2006	-	-	-	-	-
	V2007	-	-	-	+++	-