

Human MYO Mouse mAb

Cat. No: V0801
Name: Human MYO Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 8C1
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Lyophilized or Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Human Myoglobin protein, Purified native protein ([bs-41107P](#))
Pairing Information:

ELISA		Detection		
		V0801	V0802	V0803
Capture	V0801	-	+++	-
	V0802	-	-	+++
	V0803	-	-	-
ICA		Detection		
		V0801	V0802	V0803
Capture	V0801	-	+++	-
	V0802	-	-	+++
	V0803	-	-	-