

Human MYO Mouse mAb

Cat. No: V0803
 Name: Human MYO Mouse mAb
 Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
 Clone type: Monoclonal
 Clone No: 3E2
 Host/Source: Mouse
 Subtype: IgG1
 Reactivity: (predicted: Human)
 Format: Liquid
 Buffer: 0.01M PBS (pH7.4).
 Applications: ELISA, ICA
 Antigen: Human Myoglobin protein, Purified native protein ([bs-41107P](#))
 Pairing Information:

ELISA		Detection		
		V0801	V0802	V0803
Capture	V0801	-	+++	-
	V0802	-	-	+++
	V0803	-	-	-
ICA		Detection		
		V0801	V0802	V0803
Capture	V0801	-	+++	-
	V0802	-	-	+++
	V0803	-	-	-