

Human CEA Mouse mAb

Cat. No: V1906
Name: Human CEA Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 3A2
Host/Source: Mouse
Subtype: IgG2a
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ICA
Antigen: Recombinant Human CEACAM5 Protein ([bs-43141P](#))
Pairing Information:

ELISA		Detection				
		V1903	V1904	V1905	V1906	V1907
Capture	V1903	-	+++	-	-	-
	V1904	-	-	-	-	-
	V1905	-	+++	-	-	-
	V1906	-	-	-	-	-
	V1907	-	-	-	-	-
ICA		Detection				
		V1903	V1904	V1905	V1906	V1907
Capture	V1903	-	-	-	-	-
	V1904	-	-	-	-	-
	V1905	-	-	-	-	-
	V1906	-	-	-	-	-
	V1907	-	-	-	+++	-