

bs-12194R

[Primary Antibody]

HOXB6 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 3216</p> <p>Target: HOXB6</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human HOXB6: 145-224/224.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: HoxB6 is cytoplasmically expressed throughout fetal epidermal development, but displays nuclear expression in normal adult skin. HoxB6 protein also has nuclear expression in hyperproliferative skin conditions, but appears to be localized in the cytoplasm in basal and squamous cell carcinomas. HoxB6 genes are also expressed in normal adult lung.</p>	<p>Isotype: IgG</p> <p>SWISS: P17509</p> <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Cow)</p> <p>Predicted MW.: 25 kDa</p> <p>Subcellular Location: Nucleus</p>
---	---