

bs-20528R

[Primary Antibody]

CLEC2H Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 94071</p> <p>Target: CLEC2H</p> <p>Immunogen: KLH conjugated synthetic peptide derived from mouse CLEC2H: 131-218/218.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8C1T8</p>	<p>Applications: WB (1:500-2000) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Mouse)</p> <p>Predicted MW.: 24 kDa</p> <p>Subcellular Location: Cell membrane</p>
---	--	---