
CLEC2H Rabbit pAb

Catalog Number: bs-20528R

Target Protein: CLEC2H

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: (predicted:Mouse)

Predicted MW: 24 kDa

Entrez Gene: 94071

Swiss Prot: Q8C1T8

Source: KLH conjugated synthetic peptide derived from mouse CLEC2H: 131-218/218.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.