

bs-20528R

[Primary Antibody]

CLEC2H Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 94071	SWISS: Q8C1T8	Reactivity: (predicted: Mouse)
Target: CLEC2H		
Immunogen: KLH conjugated synthetic peptide derived from mouse CLEC2H: 131-218/218. < Extracellular >		
Purification: affinity purified by Protein A		Predicted MW.: 24 kDa
Concentration: 1mg/ml		Subcellular Location: Cell membrane
Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		